

# Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody Catalog # ASR1234

#### **Specification**

### Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody - Product Information

Description Anti-MOUSE IgG IgA IgM (H&L) (GOAT) (Min

X Bv, Hs, Hu Serum Proteins)

Host Goat

Conjugate Unconjugated

Target Species
Reactivity
Mouse
Application
,1,10,15,

Application Note ELISA 1:200,000; Western Blot

1:2,000-1:10,000;Immunochemistry

1:1,000-1:5,000

Physical State Liquid (sterile filtered)

Host Isotype IgG

Target Isotype IgG IgA IgM

Buffer 0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Immunogen Mouse IgG IgA and IgM whole molecule

Stabilizer None

Preservative 0.01% (w/v) Sodium Azide

# Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody - Additional Information

### **Shipping Condition**

Wet Ice

# **Purity**

Secondary Antibody was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by cross adsorption to remove unwanted specificities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgG, Mouse IgA and Mouse IgM. No reaction was observed against bovine, horse or human serum proteins. This reagent is suitable for the detection of all mouse isotypes and chain combinations. ELISA was used to confirm specificity at less than 1% of the target signal.

#### **Storage Condition**

Store vial at 4° C prior to opening. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

#### **Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

## Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody - Protein Information



# Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Anti-Mouse IgG IgA IgM (H&L) Pre-Adsorbed Secondary Antibody - Images